



Hb Hemoglobin Testing System

Accurately detect Hb and Hct levels



Fast Results



Small Sample



Portable



Sample Type:
Capillary and Venous blood

Quick

- Hemoglobin (Hb) and Hematocrit (Hct) results in < 15 seconds
- Simply insert strip, apply specimen and read results
- Low maintenance and quick clean up

Accurate

- Precise results equivalent to automatic hematology analyzers
- Wide Hb measurement range of 4.5-25.6 g/dL with Hct range of 13-75%
- Hb total precision CV $\leq 3\%$

Convenient

- Only 10 μL capillary or venous blood required
- Minimal training required, no specimen preparation
- Easy-to-read, large LCD
- Quick data transfer via Mini USB port

Specifications

FEATURE	TECHNICAL SPECIFICATION
Methodology	Reflectance Photometry
Detection Principle	Methemoglobin
Time to Results	< 15 seconds
Memory	1,000 tests with date/time and ID number
Specimen Volume	10 μL
Specimen Type	Capillary and Venous whole blood
Hb Measurement Range	4.5-25.6 g/dL (2.8-15.9 mmol/L)
Hct Range	13-75%
Wavelengths	525 nm
PC Interface	Mini USB Port
Calibration	Automatic
Hb Within Run Precision CV	$\leq 3\%$
Hb Total Precision CV	$\leq 3\%$
Accuracy	Venous Blood: $Y=0.9582X+0.5673$, $R^2=0.992$; Capillary Blood: $Y=1.006X+0.026$, $R^2=0.993$.
Operating Conditions	10-40°C (50-104°F); $\leq 90\%$ RH
Meter Storage Conditions	0-50°C (32-122°F); $\leq 90\%$ RH
Strip Storage Conditions	2-30°C (36-86°F); $\leq 85\%$ RH
Strip Shelf Life	2 years unopened canister; 3 months opened canister
Control Solution Shelf Life	6 months unopened bottle; 30 days opened bottle
Power Source	3 AAA Batteries or AC Adaptor
Battery Life	2,700 tests or 360 hours
Automatic Shut Off	8 minutes
Meter Dimensions (L X W X H)	127 mm \times 58 mm \times 25 mm (5.0" \times 2.28" \times 0.09")
LCD Dimensions (L X W)	39 mm \times 37 mm (1.54" \times 1.46")



Hb Hemoglobin Meter

Ordering Information

PRODUCT NAME	CATALOG NO.	COMPONENTS			
Hb Hemoglobin Testing System	C111-3021√	1 Meter	1 Warranty Card	1 Quick Reference Guide	2 Control Strips
		10 Lancets (26G)	10 Test Strips	1 Lancing Device Insert	1 Control Strip Insert
Hb Hemoglobin Test Strips	C131-3011√	3 AAA Batteries	1 Lancing Device	1 Test Strip Insert	10 Capillary Transfer Tubes - Plastic
		1 Users Manual	1 Carrying Case	1 Code Chip	
		25 Test Strips		1 Code Chip	1 Test Strip Insert
	C131-3021√	50 Test Strips (25/Canister)		1 Code Chip	1 Test Strip Insert
		100 Test Strips (25/Canister)		1 Code Chip	1 Test Strip Insert
Hb Hemoglobin Control Strips	C121-3031√	50 Test Strips (25/Canister)		1 Code Chip	1 Test Strip Insert
		50 Capillary Transfer Tubes - Glass/10 μL (25/Canister)			
Hb Hemoglobin Control Solution	C121-3091√	2 Control Strips (2/Canister)		1 Control Strip Insert	
		2 Bottles of Level-0 Control Solution (1 mL/bottle)			
		2 Bottles of Level-1 Control Solution (1 mL/bottle)			1 Control Solution Insert
		2 Bottles of Level-2 Control Solution (1 mL/bottle)			
Capillary Transfer Tubes	C121-3081√	50 Capillary Transfer Tubes - Plastic/10 μL			
Lancets	C121-3041*	100 Lancets (26G)(25/Bag)			
Lancing Device	C121-3051√	1 Lancing Device		1 Lancing Device Insert	
Safety Lancets I	C121-3061*	25 Safety Lancets (21G/2.8 mm); 100 Safety Lancets (28G/1.8 mm)			
Hb Hemoglobin Adaptor Kit	C121-3011√	1 Power Adaptor			

Attention: By approving the enclosed design draft, you (ACON's Customer) accept all responsibility for the accuracy of the design. If an error is detected following the printing or manufacturing of a material, you (ACON's Customer) are responsible for the cost of any inventory which is deemed unsuitable for sale and/or not compliant with local regulations .

☐ US

☒ OUS

☐ DOMESTIC

☐ OTHER

Description	Mission Hb Hemoglobin Testing System sell sheet	Part Number	1150525807	Size	210x280mm
Printing Contents		L Number		Size	
Designer	Aileen	Design Date/Version	Jul 11, 2017 /E		
Artwork checked by		Material	230g 双铜 +油性上光	Checked by	
Approved by Customer		Approved by Marketing/Sales			
Approved by P.M.T.		Approved by QA		Effective Date	